

Step One

MiniOne Gel Kit Packing Guide Crate #2

Video Tutorial




bit.ly/HowToPackMini

Parafilm and bag all liquid containers and place upright.



Stack gel boxes only **two** high!



Do not place heavy items on gel boxes.

Step Two

MiniOne Gel Kit Packing Guide Crate #2

Video Tutorial



bit.ly/HowToPackMini

Place
microfuges
on a **flat
surface** with
cords **folded
and neatly
wrapped.**

Stack
microtube
boxes **closed
and upright.**



Stack photo
hoods only
three high!



Science Education Partnership
fredhutch.org/SEP (206) 667-4487 sep@fredhutch.org



Fred Hutch
Cancer Center