

Micropipet Protocol Card - **FINNPIPETTE**

NEVER crank the micropipets above or below their ranges!!

Range:
1.0-10 μ l



Example volumes:




1.2 μ l



6.0 μ l

Range:
2-20 μ l



Example volumes:



12.5 μ l



9.75 μ l

Range:
20-200 μ l



Example volumes:



25.0 μ l



183.0 μ l

Range:
200-1000 μ l



Example volumes:



1000.0 μ l



468.0 μ l

Always hold the micropipette in a **vertical position** when there is fluid in the tip.

Tip Selection Guide

