


Range:
0.5-10 μ l



Example
 volumes:



Range:
2-20 μ l




Example
 volumes:



Range:
10-100 μ l




Example
 volumes:



Range:
20-200 μ l




Example
 volumes:



Range:
100-1000 μ l




Example
 volumes:



Range:
0.5-10 μ l




Example
volumes:



Range:
2-20 μ l




Example
volumes:



Range:
10-100 μ l




Example
volumes:



Range:
20-200 μ l




Example
volumes:



Range:
100-1000 μ l



Example
volumes:




Micropipet Volume Quiz **KEY** - Oxford

Range:
0.5-10 μ l



Example
volumes:




0.75 μ l

12.3 μ l

Range:
2-20 μ l



Example
volumes:




5.5 μ l

10.8 μ l

Range:
10-100 μ l



Example
volumes:



10.8 μ l

88.3 μ l

Range:
20-200 μ l



Example
volumes:




70.0 μ l

155.5 μ l

Range:
100-1000 μ l



Example
volumes:



100.0 μ l

1000 μ l